



Article description **Order quantity**

17.070.00 Ridge Preservation Kit acc. to Prof. Palti consisting of:		
85.185.23 Washtray for Ridge Preservation Kit		
26.690.01 Power Periostome P1 XL straight, in balanced design handle, purple color code		
26.690.02 Power Periostome P3 XL curved, mesial, in balanced design handle, blue color code		
26.690.03 Power Periostome P3 XL curved, distal, in balanced design handle, green color code		
26.182.12R Periostome P2, ZEPF-Line, rigid		
26.182.26R Periostome Combination P3 / P4, ZEPF-Line, rigid		
41.864.60 Raspatory acc. to Palti, ZEPF-Line pointed 2 mm / round & curved 2 mm		
41.864.40 Periosteal Elevator ZEPF-Line, round 4.2 mm / lance 3.8 mm	[]
46.035.05 Papillex acc. to Dr. Weng (papilla elevator), ZEPF-Line		
41.868.08M Micro Sinus Plugger, blue titanium handle, ZEPF-Line, double-ended, 17.5 cm, plugger Ø 3.5 mm and Ø 2.6 mm		
19.714.21 Augmentation Material Applicator, bucket, single-ended		
85.251.00 Mixing Cup, stainless steel, with POM lid, Ø 4 cm (kit consists of 2 pieces)	[]



Optional accessory (not included in kit) **Order quantity**

26.194.08 ZEPF BIONIK Universal Handle, PEEK, for exchangeable inserts M2.5, single-ended, green		
24.062.24 Lipcare Mouth Mirror, titanium, Ø 24 mm, Front-Surface, sterilizeable	[]
24.212.01 Combination Instrument between periodontometer North Carolina CNC (1-15 in mm-steps) & probe 3A, in ZEPF BIONIK Universal Handle, black, M4 x 0.5, exchangeable	[]



Ridge Preservation Kit



Distributed by:

Produced by:




in cooperation with Excellent Centre



videnti

As an optimal extension to this kit, you may choose between the following extraction systems

- Our X-Tool Instruments (page 05-23 to 05-25 in our catalogue)
- Our RoBa Extraction Forceps (page 05-15 to 05-17 in our catalogue)
- The Benex Extraction System (page 05-01 to 05-04 in our catalogue)



Prof. Palti

Pi
SERIES

17.070.00
Complete
Ridge Preservation Kit
acc. to **Prof. Palti**



19.714.21
Augmentation Material Applicator,
scoop, single-ended



Ridge Preservation Kit

Application



41.868.08M Micro Sinus Plugger, blue titanium handle, ZEPF-Line, double-ended,
17.5 cm, plugger Ø 3.5 mm and Ø 2.6 mm



46.035.05 Papillex acc. to Dr. Weng (papilla elevator) ZEPF-Line



41.864.40 Periosteal Elevator ZEPF-Line, round 4.2 mm / lance 3.8 mm




41.864.60 Raspatory acc. to Palti, ZEPF-Line pointed (2 mm) / round & curved (2 mm)



26.182.26R Periotome Combination P3 / P4, ZEPF-Line, rigid



26.182.12R Periotome P2, ZEPF-Line, rigid



Optional 24.212.01
Combination Instrument between periodontometer North Carolina CNC (1-15 in mm-steps) & probe 3A,
in ZEPF bIONIK Universal Handle, black, M4 x 0.5, exchangeable



26.690.01 Power Periotome P1 XL straight,
in balanced design handle, blue color code



26.690.02 Power Periotome P3 XL curved, mesial,
in balanced design handle, blue color code



26.690.02 Power Periotome P3 XL curved, distal,
in balanced design handle, green color code



reddot design award
winner 2010



Optional 26.194.08 ZEPF bIONIK Universal Handle, PEEK,
for exchangeable inserts M2.5, single-ended, green

Application



HELMUT ZEPF
MEDIZINTECHNIK GMBH

Tel.: +49 (0) 74 64 / 98 88 - 0



HELMUT ZEPF
MEDIZINTECHNIK GMBH

www.zepf-dental.com